

# 联和层析 应用文章精选

2017 | 04



联和层析致力于提供无限可能的色谱解决方案。十多年来，我们与全球多家专业色谱厂商保持良好合作伙伴关系，并将他们优势产品引进到国内。本着**诚恳、实在、专业、负责**的理念，联合层析在样品前处理及自动化的领域已占据领先地位。



为服务广大客户对于了解国际间最新色谱应用方向的需求，联和层析定期精选并翻译应用文章，为国内的科技服务提升，略尽绵薄之力！

## 《透过 PAL SPME Arrow 和气相色谱/质谱仪自动化测定水中的亚硝胺》

Quantification of nitrosamines in water by automated PAL SPME Arrow and GC/MS ; 原文刊载于 CTC 于 2016 年 2 月发行的 GC Application Note ; 由联和层析贸易 (上海) 有限公司编译

### 【应用背景介绍】

对于许多动物 (包含人类) 而言, 大多数亚硝基化合物是可疑致癌物。人类可能透过摄入或吸入环境中的化合物而接触致癌的亚硝胺。这些亚硝胺是由与胺类化合物反应形成, 尤其是二级胺和三级胺或含亚硝酸盐的胺基化合物 (例如: 食物中的防腐剂)。NDMA 和 NDEA 是食品材料中最常见的亚硝胺。



## 与我们联系



联和层析贸易 (上海) 有限公司

www.unichrom.cn

上海市浦东新区向城路 19 号 10A

电话: 86-21-64479550

传真: 86-21-20227695

邮箱: info@unichrom.cn

北京市朝阳区锦芳路 1 号 9 号楼 609 室

电话: 86-10-84869268

传真: 86-10-84869268

邮箱: james-jiang@unichrom.cn

广州市天河区天坤三路和美街 11 号 1706 室

电话: 15217202438

邮箱: gz@unichrom.cn

### 【结论】

- ◆ 制定出一个为多种亚硝胺于水溶液样品中的定量分析方法
- ◆ 根据先前的研究, 选定 Carboxen / PDMS SPME 纤维为吸附物
- ◆ Carbon WR SPME Arrow 提取率/检测器信号比 Carboxen / PDMS SPME 纤维提升 3-5 倍
- ◆ 检测限 ( =S/N>3 ) 范围从 5ng/L ( N-nitrosodipropylamin ) 到 50N/L ( n-nitrosopyrolindinamine )
- ◆ 在 1µg/L 的重复性为 4-12%
- ◆ 运行一个样品的时间为 30 分钟提取+30 分钟 GC 运行时间。重叠提取 7 个样本的总运行时间约为 240 分钟
- ◆ PAL RTC/RSI 手臂配置 SPME Arrow 模块, 可全自动进行 SPME 分析, 并确保过程的安全及结果的高重复性

SPME Arrow

萃取

解析/注射

调节

# GC Application Note

Quantification of nitrosamines in water by automated PAL SPME-Arrow and GC/MS



FOOD SAFETY





# Quantification of nitrosamines in water by automated PAL SPME-Arrow and GC/MS


Peter Egli, Beat Schilling, BGB Analytik AG, Adliswil, Switzerland  
Guenter Boehm, CTC Analytics AG, Zwingen, Switzerland

## Short summary:

- A method for the quantitative analysis of a number of different N-nitrosamines from aqueous samples has been worked out.
- Carboxen/PDMS SPME fibers and Carbon WR SPME Arrow were compared. Carbon WR SPME Arrow gave a 3-5 x higher extraction yield/detector signal.
- The limits of detection (= S/N > 3) range from 5 ng/L for N-nitroso-dipropylamin to 50 n/L for N-nitroso-pyrolindin-amine.
- Repeatability @ 1 µg/L ranges from 4-12%.

## Introduction:

Most nitroso compounds are suspected carcinogens in various animals including human beings. Human exposure to carcinogenic N-nitrosamines may result directly from ingestion or inhalation of preformed compounds from the environment. They are formed by the reaction of amines, especially secondary and tertiary amines or amino group containing compounds with nitrite (present e.g. as food preservative). N-nitroso-dimethylamine (NDMA) and N-nitroso-diethylamine (NDEA) are the most common nitrosamines found in food materials.



## Experimental:

### 1. Immersion Extraction

Sample	10ml water in 20ml HS vial with 4g KCl
Fiber	Carbon WR SPME Arrow (PAL Sytem No: ARR11-C-WR-120-20-P1) Carboxen/PDMS SPME Fiber ( part of Supelco kit, No: 57284-U)
Pre conditioning	0:30 min
Pre incubation time	1:00 min
Incubation temp	50°C
Agitation speed	500 rpm
Needle penetration	22 mm
Fiber penetration	30 mm
Extraction time	30 min
Desorption time/temp	2:00 min/ 260°C

### 2. Chemicals

Nitrosamines calibration standard (Restek No: 31898)

N-nitroso-dimethylamine

N-nitroso-methylethylamine

N-nitroso-diethylamine

N-nitroso-dipropylamine

N-nitroso-dibutylamine

N-nitroso-pyrolidone

N-nitroso-piperidine

Diluted in Methyl tert-butyl ether (Fluka No: 20249)

### 3. GC/MS

GC	Varian 3400
MS	Varian Saturn ion trap
Column	60 m; 0.25 mm ; 25 µm BGB Wax 20M
Carriergas	Hydrogen 10.0 psi
Temp. program	50°C for 1min, then 10°C/min to 240°C
Injector	260°C
Mass Range	42 bis 250 m/z

With the PAL RTC and RSI the entire SPME process is fully automated. This guarantees process safety and high reproducibility.



**Overlaid Chromatogram Plots**

Plot 1: c:\... \desktop\spme\nitrosamine\iam in 22 .sm s Ions: 74+88+102+102+130+84+114+100 all

Plot 2: c:\... \desktop\spme\nitrosamine\iam in 26 .sm s Ions: 74+88+102+102+130+84+114+100 all




Fig. 1: Chromatograms of 1 µL of a 1 µg/L standard solution and blank direct injection





Fig. 2: Comparison of relative detector (MS) response for different nitrosamines @ 1 µg/L between Carbon WR SPME Arrow and Carboxen/PDMS SPME

	N-nitroso-dimethylamine	N-nitroso-methylethylamine	N-nitroso-diethylamine	N-nitroso-dipropylamine	N-nitroso-dibutylamine	N-nitroso-piperidine	N-nitroso-pyrrolidine
Average (n=5)	153621	1105607	430555	566316	218648	218468	264031
% RSD	11.9	8.8	5.9	4.2	6.8	14.1	10.4


Tab. 1: Reproducibility for the quantification of different nitrosamines @ 1 µg/L with Carbon WR SPME Arrow.

Overlaid Chromatogram Plots




N-nitroso-diethylamine @ 50 ng/L, S/N 24

Overlaid Chromatogram Plots



N-nitroso-dipropylamine @ 50 ng/L, S/N 34

Overlaid Chromatogram Plots



N-nitroso-pyrrolidinamine @ 50 ng/L, S/N 3 (=LOD)

Fig. 3: Chromatograms for selected N-nitrosamines, @ 1000 ng/L and 50 ng/L with Carbon WR SPME Arrow.






Fig. 4, 5, 6: Calibration curves for N-nitroso-diethylamine (top), N-nitroso-dipropylamine (middle) and N-nitroso-pyrrolidineamine (bottom)

## Conclusions:

- A method for the quantitative analysis of a number of different N-nitrosoamines from aqueous samples has been worked out.
- Based on previous publications (ref.1) Carboxen/PDMS SPME fibers have been selected as sorbent material
- Carboxen/PDMS SPME fibers and Carbon WR SPME Arrow were compared. Carbon WR SPME Arrow gave a 3-5 x higher extraction yield/detector signal.
- The limits of detection (= S/N > 3) range from 5 ng/L for N-nitroso-dipropylamine to 50 ng/L for N-nitroso-pyrrolidineamine.
- Repeatability @ 1 µg/L ranges from 4-12% @ 1 µg/L.
- Time for one sample was 30 min extraction + 30 min GC runtime. With overlapped extraction (e.g. with PAL Sample Control <http://www.palsystem.com/index.php?id=243>, or many other chromatographic data systems) the total runtime for 7 samples was approximately 240 min.

## References:

- [1] Grebel JE, Young CC, Suffet IH, J. Chrom. A, 1117 (2006) 11-18

---

## Legal Statements

CTC Analytics AG reserves the right to make improvements and/or changes to the product(s) described in this document at any time without prior notice.

CTC Analytics AG makes no warranty of any kind pertaining to this product, including but not limited to implied warranties of merchantability and suitability for a particular purpose.

Under no circumstances shall CTC Analytics AG be held liable for any coincidental damage or damages arising as a consequence of or from the use of this document.

© 2016 CTC Analytics AG. All rights reserved. Neither this publication nor any part hereof may be copied, photocopied, reproduced, translated, distributed or reduced to electronic medium or machine readable form without the prior written permission from CTC Analytics AG, except as permitted under copyright laws.

CTC Analytics AG acknowledges all trade names and trademarks used as the property of their respective owners.

PAL is a registered trademark of CTC Analytics AG | Switzerland

## Imprint

Date of print: 02.2016

CTC Analytics AG  
Industriestrasse 20  
CH-4222 Zwingen  
Switzerland  
T +41 61 765 81 00  
Contact: [info@ctc.ch](mailto:info@ctc.ch)